

VICTO

Adjustable Artificial Sphincter

ORDERING INFORMATION

Starter Kits

DESCRIPTION	CUFF SIZE	CODE
VICTO STARTER KIT 3.7	3.7	KITVICTO-37
VICTO STARTER KIT 4.0	4.0	KITVICTO-40
VICTO STARTER KIT 4.5	4.5	KITVICTO-45
VICTO STARTER KIT 5.0	5.0	KITVICTO-50
VICTO STARTER KIT 5.5	5.5	KITVICTO-55

Accessories Kit

DESCRIPTION	CODE
VICTO ACCESSORIES KIT	KITAC-02

victosphincter.com

Magnetic Resonance (MR) Environment

Non-clinical testing has demonstrated that Victo is MR conditional.

A patient with this device can be scanned safely under the following conditions:



- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of < 2 W/kg (Normal Operating Mode)

Promedon
Urology and Pelvic Floor

Company information: Distribution EMEA & APAC | Promedon GmbH | orders@promedon.com | Phone: +49 8031 900 400 | Fax: +49 08031 900 4040 | An der Alten Spinnerei 5 | 83059 Kolbermoor | GERMANY

Manufactured by: PROMEDON S.A. | www.promedon-upf.com



B-31-11 / PSD-00000317/011 / MAR-2022

VICTO

Adjustable Artificial Sphincter



TAKING THE ARTIFICIAL URINARY SPHINCTER CONCEPT TO THE NEXT LEVEL


- Minimally invasive adjustment procedure
- Pre-connected device
- Ergonomic pump design

VICTO ARTIFICIAL URINARY SPHINCTER.

VICTO is a preconnected artificial urinary sphincter that enables pressure adjustments with a minimally invasive ambulant procedure.

MINIMALLY INVASIVE ADJUSTMENT PROCEDURE

With its unique pump design, the system pressure in VICTO can easily be adjusted through the self-sealing port, in an outpatient and easy-to-perform procedure.



ERGONOMIC PUMP DESIGN

The soft silicone pump has an easy-to-palpate design, which makes it simple for the patient to control. With only one pump stroke the urethra is decompressed and micturition is enabled.



PRE-CONNECTED DEVICE

VICTO is provided pre-connected, ready for filling and implantation.



VICTO is available in 5 different cuff sizes, which can be selected intraoperatively using the cuff size selection tape, provided in the Accessories kit.

3.7_{CM}

4.0_{CM}

4.5_{CM}

5.0_{CM}


5.5_{CM}

ACTIVATION & ADJUSTMENT PROCEDURE

Preparation during surgery: +13ml

Activation +4ml

Individual adjustments(+1 ml), if needed



NEEDLE RECOMMENDATION:

The usage of a non-coring needle 25G is recommended for the filling or liquid removal procedures performed during the Victo activation, adjustment or deactivation.

INTENDED USE

VICTO is intended for the reestablishment of urinary continence in the defined patient population.

INDICATIONS

VICTO is indicated for the treatment of male stress urinary incontinence.

CONTRAINDICATIONS

- Uncontrolled overactive bladder
- Untreated or irreversible infravesical obstruction
- Impaired cognitive or manual abilities which may affect manual dexterity or motivation, and prevent the patient from operating the device
- Untreated acute urinary tract infection or active infection at the surgical sites or systemic infections (sepsis)
- Known allergy to silicone or titanium
- Patients with stricture/stenotic disease
- Patients whom the surgeon determines to be not suitable due to risks associated with open surgical procedures, and/or with the patients' medical history (physical or mental problems)
- Prior abdominal surgery, status post peritonitis, and other conditions resulting in peritoneal adhesions, fibrosis, scars that could potentially interfere with the correct functioning of the pressure regulating balloon

ORDERING INFORMATION

Starter Kits		
DESCRIPTION	CUFF SIZE	CODE
VICTO STARTER KIT 3.7	3.7	KITVICTO-37
VICTO STARTER KIT 4.0	4.0	KITVICTO-40
VICTO STARTER KIT 4.5	4.5	KITVICTO-45
VICTO STARTER KIT 5.0	5.0	KITVICTO-50
VICTO STARTER KIT 5.5	5.5	KITVICTO-55

Accessories Kit	
DESCRIPTION	CODE
VICTO ACCESSORIES KIT	KITAC-02

victosphincter.com

Magnetic Resonance (MR) Environment

Non-clinical testing has demonstrated that Victo is MR conditional. A patient with this device can be scanned safely under the following conditions:



- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of < 2 W/kg (Normal Operating Mode)



TAKING THE ARTIFICIAL URINARY SPHINCTER CONCEPT TO THE NEXT LEVEL

- Minimally invasive adjustment procedure
- Pre-connected device
- Ergonomic pump design

ORDERING INFORMATION

Starter Kits

DESCRIPTION	CUFF SIZE	CODE
VICTO STARTER KIT 3.7	3.7	KITVICTO-37
VICTO STARTER KIT 4.0	4.0	KITVICTO-40
VICTO STARTER KIT 4.5	4.5	KITVICTO-45
VICTO STARTER KIT 5.0	5.0	KITVICTO-50
VICTO STARTER KIT 5.5	5.5	KITVICTO-55

Accessories Kit

DESCRIPTION	CODE
VICTO ACCESSORIES KIT	KITAC-02

victosphincter.com

Magnetic Resonance (MR) Environment

Non-clinical testing has demonstrated that Victo is MR conditional. A patient with this device can be scanned safely under the following conditions:



- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of < 2 W/kg (Normal Operating Mode)

Promedon
Urology and Pelvic Floor



TAKING THE ARTIFICIAL URINARY SPHINCTER CONCEPT TO THE NEXT LEVEL

- Minimally invasive adjustment procedure
- Pre-connected device
- Ergonomic pump design

ORDERING INFORMATION

Starter Kits

DESCRIPTION	CUFF SIZE	CODE
VICTO STARTER KIT 3.7	3.7	KITVICTO-37
VICTO STARTER KIT 4.0	4.0	KITVICTO-40
VICTO STARTER KIT 4.5	4.5	KITVICTO-45
VICTO STARTER KIT 5.0	5.0	KITVICTO-50
VICTO STARTER KIT 5.5	5.5	KITVICTO-55

Accessories Kit

DESCRIPTION	CODE
VICTO ACCESSORIES KIT	KITAC-02

victosphincter.com

Magnetic Resonance (MR) Environment

Non-clinical testing has demonstrated that Victo is MR conditional. A patient with this device can be scanned safely under the following conditions:



- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of < 2 W/kg (Normal Operating Mode)

Promedon
Urology and Pelvic Floor



TAKING THE ARTIFICIAL URINARY SPHINCTER CONCEPT TO THE NEXT LEVEL

- Minimally invasive adjustment procedure
- Pre-connected device
- Ergonomic pump design