inDrive

70-15-3219 (02) / 23-Jun-2021

Instructions for use
ENGLISH

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Promedon

INSTRUCTIONS FOR USE

DESCRIPTION

The InDrive ureteral access sheath is made up of two components:

1. The sheath itself with a connector.
2. Dilator with a connector.

PACKAGE CONTENTS

InDrive is made up of both parts mentioned in the description (dimensions indicated on the label).

INDICATIONS

This device is used to maintain a working channel open during endoscopic urological procedures with the aim of facilitating the passage of the endoscope and other instruments into the urinary tract.

CONTRAINDICATIONS

Impossibility of or excessive strength in placing the sheath.

PRECAUTIONS

Although the ureteral access sheath may be inserted with no additional maneuvering, difficulty may be encountered in its passage. Certain groups have been identified as at risk:

- Young males with narrow ureters.
- Duplicate collecting systems.
- Small caliber ureters.

Single-use sterile product (sterilized with ethylene oxide), do not reuse or resterilize, since this could reduce the performance of the device and increase the risk of improper sterilization and cross-contamination.

Do not use kit if the package is open or damaged.

Only physicians familiar with the ureteral endoscopic procedure are qualified to use it.

STORAGE

It is advisable to store the kit in its original packaging at: Room TEMPERATURE.

USE BEFORE THE EXPIRY DATE SPECIFIED ON THE PACKAGING.


POSSIBLE COMPLICATIONS

- Perforation of the ureter.

PROCEDURE

- Ensure the sheath-dilator set is properly attached.
- Activate the hydrophilic coating, washing the sheath-dilator set with isotonic saline.
- Once the guide (preferably 0.035" in diameter) has been inserted, move the device tip forward (where the dilator is conical in shape) on the guide to the interior of the ureter.
- Use fluoroscopy to confirm that both the dilator and the sheath are correctly placed.
- Unlock the connector of the dilator from the connector of the sheath and remove the dilator and guide.

SYMBOLS USED IN LABELS

| | |
|---|--|
|  | CATALOG NUMBER |
|  | LOT NUMBER |
|  | CAUTION |
|  | EXPIRY DATE |
|  | DO NOT REUSE |
|  | DATE OF MANUFACTURE |
|  | STERILE. STERILIZATION METHOD: ETHYLENE OXIDE |
|  | DO NOT RESTERILIZE |
|  | MANUFACTURER |
|  | SERIAL NUMBER |
|  | CONSULT INSTRUCTIONS FOR USE |
|  | DO NOT USE IF PACKAGING IS DAMAGED |
|  | KEEP AWAY FROM SUNLIGHT |
|  | KEEP DRY |

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Fabricado por: PROMEDON S.A
PROMEDON S.A.