PERKUT

DILATOR SET FOR PERCUTANEOUS
NEPHROLITHOTOMY

M-70-15-3607 (02) / 27 - May - 2021

Instructions for use
ENGLISH

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Promedon

INSTRUCTIONS FOR USE

PACKAGE CONTENTS AND DESCRIPTION

The set of dilators for percutaneous nephrolithotomy **PERKUT (ST-PCNL-LB)** includes:

- 13 dilators (6 Fr–30 Fr)
- 1 introducer catheter (8 Fr)
- 1 introducer catheter sheath (10 Fr)
- 4 sheaths (24 Fr–30 Fr internal diameter)

INDICATIONS FOR USE

Dilate percutaneous tract before performing the percutaneous renal endoscopy in order to introduce the nephroscope and its accessories.

CONTRAINDICATIONS

This product is contraindicated in the presence of conditions which create unacceptable risks during catheterization.

Previous kidney surgery.

PRECAUTIONS

Single-use sterile product (sterilized with ethylene oxide). Any institution, professional, or third party who reprocesses, reconditions, remanufactures, sterilizes and/or reuses these disposable devices shall be entirely responsible for their safety and efficacy.

Promedon has not designed these devices to be reprocessed or reused. Therefore, cannot verify that the reprocessing ensures the cleaning and/or sterilization of any device, or the maintenance of its structural integrity, in order to ensure patient and/or user safety.

Operating room conditions must be in accordance with hospital, administrative, or local government procedures.

After use, dispose of the product and packaging in accordance with hospital, administrative, or local government procedures.

Do not use device if the package (sterile barrier) is opened or damaged.

This device should not be in contact with any electrical instrument.

Only a physician skilled in percutaneous renal endoscopic procedures is competent in the proper use of this device.

STORAGE

It is advisable to store the device in its original packaging at ROOM TEMPERATURE.

USE BEFORE THE EXPIRY DATE ON THE PACKAGING.

POTENTIAL COMPLICATIONS

- Ureteral orifice injury
- Hematuria
- Delayed stenosis

Adverse events that might result from a dilation procedure include:

- Tissue trauma
- Tissue perforation

PROCEDURE

PREPARATION:

Start sequential dilation using 6 Fr to 10 Fr dilators after inserting the guidewire in the percutaneous puncture tract.


Insert the introducer catheter into the guidewire. Then, position the introducer catheter sheath in the introducer catheter.

Complete system dilation up to the required size (24 Fr, 26 Fr, 28 Fr or 30 Fr).

Slide the corresponding sheath on the dilator to create the working channel.

Remove dilator and leave sheath in place.

SYMBOLS USED IN LABELS

	CATALOG NUMBER
	LOT NUMBER
	SERIAL NUMBER
	CAUTION
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	READ THE INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	DO NOT RESTERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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Fabricado por: PROMEDON S.A