URODRAIN

PSD-00000354 (03) / AUG-2022

Instructions for use
ENGLISH

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Promedon

Promedon

DESCRIPTION

Urodrain is a single-use product for temporary drainage between the pyeloureteral junction and the bladder.

PACKAGE CONTENTS

The kit contains:

- 1 radiopaque polyurethane stent with/without hydrophilic sheath, double pigtail and thread.
- 1 radiolucent stent pusher.
- 1 clamp (only for closed proximal end model).

There are two versions that, besides including the previously mentioned elements, contain:

- 1 PTFE coated stainless steel guide wire or
- 1 Hydrophilic Nitinol guide wire.

INDICATIONS

Urodrain is indicated for temporary drainage between the pyeloureteral junction and the bladder.

CONTRAINDICATIONS

- Urodrain must not be used in patients with severe coagulopathies.
- Urodrain may become obstructed in the presence of mycotic or fungal infections.

WARNINGS

The implant must not be handled with pointed, serrated, or sharp objects since any collapse, damage, perforation, or tearing can cause subsequent complications. Fluff, fingerprints, talcum powder, or other contaminated surfaces may cause reactions to foreign bodies. Maximum precautions must be taken to avoid contamination.

Any bending or twisting of the stent, prior or during insertion, may damage its integrity.

The use of this product in pregnant women must be monitored more frequently, because calcium incrustations may build up.

PRECAUTIONS

This product has not been designed for use as a permanent device. Cystoscopy or radiography examinations are recommended for periodic monitoring.

Urodrain components are designed for SINGLE use only. Therefore, DO NOT REUSE or RESTERILIZE, since this could reduce the performance of the device and increase the risk of improper sterilization and cross-contamination. Urodrain is provided sterile and pyrogen-free. Its packaging is a double layered protective pouch. IF ANY OF THESE ARE DAMAGED, DO NOT IMPLANT THE STENT.

Only physicians familiar with this application are qualified to use it.

STORAGE

Storage should meet the following conditions:

- TEMPERATURE: Room temperature

USE BEFORE THE EXPIRY DATE SPECIFIED ON THE PACKAGING.

COMPLICATIONS

Possible complications associated with the use of a stent should be discussed with the patient before surgery.

The use of a stent may cause associated complications related to the medication and methods used in the surgical procedure, as well as complications related to each patient's specific degree of intolerance to the implantation of any foreign body. Some complications may require stent removal.

Potential complications associated to ureteral stents are:

- Ureteral reflux.
- Extravasation.
- Catheter occlusion.
- Catheter displacement.
- Bleeding.
- Sepsis.
- Perforation of kidney, renal pelvis, ureter, and bladder.
- Peritonitis.
- Hematuria.
- Pain in the lower abdomen.

Promedon requires medical professionals to report any complication associated with the use of the ureteral stent to the Company.

PROCEDURE

Open proximal end ureteral stent:

1. Wash the ureteral stent, the pusher, and the guide with physiological saline solution.
2. Identify the flexible end of the guide.
3. Insert the flexible end of the guide through the working channel of the endoscope to be used.
4. Locate the ureter and insert the guide until the proximal end is placed in the renal pelvis.
5. Insert the ureteral stent concentrically into the previously inserted guide.
6. Push the ureteral stent in the direction of the renal pelvis with the aid of the pusher and the guide already in place, until the proximal end of the ureteral stent sets itself free, and regains its original bend.
7. Hold the pusher in this position, and withdraw the guide until the distal pigtail of the ureteral stent sets itself free inside the bladder.


Closed proximal end ureteral stent:

1. Wash the ureteral stent, the pusher, and the guide with physiological saline solution.
2. Identify the rigid end of the guide.
3. Insert the rigid end of the guide inside the ureteral stent until the pigtail becomes straight.
4. Insert the pusher until it reaches the ureteral stent.
5. Secure the pusher with the help of a clamp.
6. Push the guide mounted ureteral stent, until it reaches the renal pelvis through the working channel of the endoscope in use.
7. Opening the clamp must be done only after checking its correct positioning using X-rays.
8. Hold the pusher in this position and withdraw the guide until the distal pigtail sets itself free inside the renal pelvis, and the distal stent of the ureteral stent sets itself free inside the bladder.

Note: The ureteral stent's position may be corrected using the thread from the distal end of the pigtail, both for open or closed end ureteral stents. After the stent has been correctly positioned, cut and withdraw knot in thread.

In order achieve the hydrophilic effect, stents with sheaths must be immersed in saline solution before use.

SYMBOLS USED IN LABELS

| | |
|---|---|
|  | CATALOG NUMBER |
|  | LOT NUMBER |
|  | SERIAL NUMBER |
|  | CAUTION |
|  | EXPIRY DATE |
|  | DO NOT REUSE |
|  | DATE OF MANUFACTURE |
|  | STERILE. STERILIZATION METHOD: ETHYLENE OXIDE |
|  | MANUFACTURER |
|  | CONSULT INSTRUCTIONS FOR USE |
|  | DO NOT USE IF PACKAGING IS DAMAGED |
|  | DO NOT RESTERILIZE |
|  | KEEP AWAY FROM SUNLIGHT |
|  | KEEP DRY |

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Fabricado por: PROMEDON S.A
PROMEDON S.A.