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## INSTRUCTIONS FOR USE

### DESCRIPTION

The Promedon Nephrostomy Kit includes:

- 1 catheter
- 3 dilators
- 1 IP needle
- 1 Chiba needle
- 1 guidewire
- 1 connection adapter
- 1 stylet

### INDICATIONS OF USE

The kit is used for placing a catheter percutaneously in the renal pelvis for nephrostomy drainage.

### CONTRAINDICATIONS

They have not been described.

### WARNINGS

The catheter should not be handled with pointed, serrated or sharp objects, since any depression, deterioration, puncture or tear may cause subsequent complications. Fluff, fingerprints, powder or any other contaminated area may cause reactions to foreign bodies. Every precaution must be taken to avoid contamination. The bending or torsion of the catheter before or during placement can damage its integrity.

### PRECAUTIONS

- Promedon did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety
- Single-use sterile product (sterilized with ethylene oxide). Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, sterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.
- Operating room conditions must be in accordance with local hospital, administrative or government procedure.
- After use, dispose of the product and packaging. It is recommended to perform in accordance with local hospital, administrative or government procedure.
- Do not use the device if the packaging (sterile barrier) is open or damaged.
- This device must not be in contact with any electric instruments.
- Only a physician familiar with the nephrostomy procedure is qualified to use this device.
- Periodic evaluations are recommended. The catheter must not remain implanted for more than one month.

### STORAGE

Store the device in its original packaging at ROOM TEMPERATURE.  
DO NOT USE AFTER THE EXPIRATION DATE ON THE PACKAGING.

### PROCEDURE

1. First, locate the kidney to be drained by means of regular X-ray, intravenous pyelography, retrograde pyelography, ultrasound or CT scan. Using fluoroscopic guidance, identify the area of skin overlying the collecting system and inject a local anesthetic to numb the area.
2. Insert one of the needles into the collecting system at a right angle to the skin. A urine sample can be collected, if desired. Visualization of the pyelocaliceal system can be enhanced by injecting a contrast medium.
3. Pass one of the needles vertically through a small skin incision (made with the stylet or scalpel) and insert it in the proper section of the pyelocaliceal system. Upon removal of the obturator, the free flow of liquid through the needle confirms the satisfactory location of the tip of the needle's cannula.

4. Pass the flexible end of the guidewire through the connector of the needle cannula and move it forward into the collecting system. Confirm the position under fluoroscopic guidance. During removal of the needle's cannula, maintain the position of the guidewire over it.
5. In order to facilitate the insertion of the nephrostomy catheter, dilate the musculofascial tract using the dilators, which are provided sequentially (in increasing order), while advancing the dilator over the guidewire. Before inserting the catheter, it must be unrolled.
6. Pass the end of the catheter over the external end of the guidewire, and gradually insert the end deep into the collecting system. Confirm the position under fluoroscopic guidance. While using one hand to firmly hold the catheter's body in position, use the other hand to pull out the guidewire.
7. A dressing may be applied. Use the connecting tube to connect the catheter to a drainage bag or a leg (urine collection) bag.

### SYMBOLS USED IN LABELS

	PRODUCT CODE
	LOT NUMBER
	SERIAL NUMBER
	CAUTION
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	READ THE INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	DO NOT RESTERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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**PROMEDON S.A.**