FLUENZ

Hybrid Guidewire

M-70-15-3639 (01) / 28-May-2020

Instructions for use
ENGLISH

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Promedon

INSTRUCTIONS FOR USE

DESCRIPTION AND PACKAGING CONTENTS

The **FLUENZ (G-35-150S-COMBI)** hybrid guidewire package contains:
1 Hybrid Guidewire.

INDICATIONS FOR USE

The guidewires are intended for use in facilitating the placement of endourological devices during diagnostic or interventional procedures. These guidewires are not intended for coronary arteries, vascular or neurological use.

CONTRAINDICATIONS

These guidewires are not intended for coronary arteries, vascular or neurological use.

For specific additional contraindications, check the label of the device that will be used with the guide.

WARNINGS

- A thorough understanding of the technical principles, clinical applications, and risks associated with the use of guidewires is necessary before using this product.
- This device is restricted to use by or under the supervision of physicians trained in urological endoscopic procedures. Exercise caution to prevent perforation or trauma of the linings and associated tissues, channels, and ducts.
- Failure to respect the previously mentioned precautions may result in injuries to the channels or ducts, abrasion of the hydrophilic coating, release of plastic fragments from the guidewire, or damage to or separation of the guidewire, which may require intervention.
- Do not withdraw the guidewire through a metal cannula or needle. Removal through a metal device may lead to the release of wire fragments in the urinary system and damage to or separation of the outer polymer jacket, which would require extraction. Extreme caution should be exercised when using a single-wall puncture needle.
- Use extreme caution when using a laser or electrocautery device in order to avoid contact with the guidewire. Direct contact may damage or sever the wire.
- Do not reshape the guidewire by any means. Attempting to reshape the wire may lead to damage and release of wire fragments in the urinary system.
- When exchanging or removing a catheter using the guidewire, secure and maintain the guidewire in place under fluoroscopy to avoid unexpected guidewire advancement. Otherwise the wire's tip may cause injuries to the urinary tract.
- Move the guidewire slowly and carefully in the urinary system while confirming the behavior and location of the wire's tip under fluoroscopy. Excessive handling of the guidewire without fluoroscopic confirmation may result in perforation or trauma of the linings or associated tissues, channels, or ducts. If there is any resistance or if the tip's behavior or location seem inadequate, STOP moving the guidewire or catheter and determine the cause by fluoroscopy. Failure to exercise proper caution may result in bending, kinking, separation of the guidewire's tip, damage to the catheter, or injuries to the urinary system. If necessary, remove the guidewire and auxiliary device or scope as a complete unit to avoid complications.
- Do not attempt to use the guidewire if it has been bent, kinked, or damaged. Use of a damaged wire may result in damage to the linings and associated tissues, channels, or ducts, or release of wire fragments in the urinary system.

PRECAUTIONS

- Single-use sterile product (sterilized with ethylene oxide). Any institution, professional, or third party that reprocesses, reconditions, remanufactures, sterilizes and/or reuses these disposable devices shall be entirely responsible for their safety and efficacy.
- Promedon has not designed these devices to be reprocessed or reused. Therefore, they cannot verify that the reprocessing ensures the cleaning and/or sterilization of any device, or the maintenance of its structural integrity, in order to ensure patient and/or user safety.
- Operating room conditions must be in accordance with hospital, administrative, or local government procedures.

- After use, dispose of the product and packaging. Do so in accordance with hospital, administrative, or local government procedures.
- Do not use the device if packaging (sterile barrier) is open or damaged.
- This device should not be in contact with any electric instrument.
- Only a physician skilled in endoscopic kidney stone removal procedures is competent to use it.
- Do not use this product without having previously read and understood all of the instructions included here.
- The entire operation should be performed in a sterile field.
- The product will be sterile in an unopened unit packaging that has not been damaged. Do not use if the unit packaging or the guidewire are broken, damaged, or soiled. The products should be used immediately after opening the package. All guidewires should be disposed of safely and properly after use, following local regulations for medical waste management.
- When using a medication or device concurrently with FLUENZ hybrid guidewires, or any other guidewire, the operator should have a full understanding of the characteristics/properties of such medication or device to avoid damaging the wire.
- The guidewire's surface is not slippery unless it is wet. Before using and inserting the wire through a catheter, fill the dispensing container and the catheter with physiological saline solution.
- The guidewire should be advanced through the scope using short movements to prevent unintended damage to the device or injury to the patient.
- When reinserting the guidewire back into the dispensing container, be careful not to damage the wire's coating with the edge of the dispensing container.
- Do not use a metal torque device with the guidewire. Using a metal torque device may damage the wire. Do not slip a tightened torque device over the wire either, as this may damage the wire.
- Due to variations of certain catheter tip inner diameters, abrasion of the coating may occur during handling. If there is any resistance while introducing the catheter, it is advisable to stop using those catheters.
- After removing it from the patient's urinary system and before reinserting it in the same patient during the same catheterization, the guidewire must be rinsed in a container filled with physiological saline solution. Use of alcohol or antiseptic solutions or other solvents must be avoided, as they may adversely affect the surface of the guidewire.

STORAGE

Store the device in its original packaging at ROOM TEMPERATURE. DO NOT USE AFTER THE EXPIRATION DATE ON THE PACKAGING.

POTENTIAL COMPLICATIONS

Complications that may arise from the use of guidewires in urological procedures include:

- Urinary tract perforation
- Acute bleeding
- Hemorrhage
- Tissue trauma
- Edema
- Foreign objects in the body
- Infection
- Hemoglobinuria
- Peritonitis
- Ureteral avulsion

PROCEDURE

PREPARATION:

Carefully examine the hybrid guidewire pouch to check for any damages that may have occurred during shipment. Hybrid guidewires are delicate instruments and should be handled with caution. Prior to use, and whenever possible during the procedure, carefully inspect the hybrid guidewire to check for bends or kinks.

Before removing the guidewire from the dispensing container (by its distal end), inject physiological saline solution into the Luer-Lock connector of the dispensing container in order to completely wet the surface of the hybrid guidewire. Use of alcohol or antiseptic solutions

or other solvents must be avoided, as they may adversely affect the surface of the hybrid guidewire. Do not use a dry cloth on the hybrid guidewire. Before inserting the guidewire through a catheter, fill the catheter with physiological saline solution.

RECOMMENDED TECHNIQUE

1. Remove the guidewire from the dispensing container. Save the container to store the hybrid guidewire if you are going to use it again during the same intervention on the same patient.
2. Prior to use, inspect the hybrid guidewire to check for the following:
 - a. Roughness or abrasion of the tip.
 - b. Kinks throughout the length of the guidewire.
3. Insert the hybrid guidewire through the natural or developed passage according to the selected surgical procedure.
4. Confirm placement of the hybrid guidewire using the selected imaging system.
5. Keep the hybrid guidewire in position while performing the procedure.
6. Slowly and carefully remove the hybrid guidewire from the patient to prevent kinks.

SYMBOLS USED IN LABELS

PRODUCT CODE

BATCH NUMBER

SERIAL NUMBER

PRECAUTION

EXPIRATION DATE

DO NOT REUSE

DATE OF MANUFACTURE

STERILE. STERILIZATION METHOD: ETHYLENE OXIDE

MANUFACTURER

CHECK INSTRUCTIONS FOR USE

DO NOT USE IF THE PACKAGING IS DAMAGED

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PROMEDON S.A.