STRATONE

FRONTAL

M-70 -15-3500/ PSD-00000258 (03) / FEB-2022

Instructions for use
ENGLISH

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Promedon

INSTRUCTIONS FOR USE**DESCRIPTION**

The universal Nitinol stone retrieval basket is comprised of three main components:

1. Stone retrieval basket.
2. Sheath.
3. Detachable handle.

PACKAGING CONTENTS

Stratone Frontal is comprised of:

- 1 (one) stone retrieval basket.
- 1 (one) Luer-lock valve.

INDICATIONS FOR USE

Stratone Frontal are indicated for the endoscopic removal of urological stones.

CONTRAINDICATIONS

- They have not been described.

PRECAUTIONS

- Some urological stones may be too large for this type of procedure.
- Promedon did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- Sterile product (sterilized with ethylene oxide), single use. Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, sterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.
- Operating room conditions must be in accordance with local hospital, administrative or government procedure.
- After use, it is recommended to discard the product and packaging. It is recommended to perform in accordance with local hospital, administrative or government procedure.
- Do not use the device if the packaging (sterile barrier) is open or damaged.
- Only a physician familiar with the endoscopic urological stone elimination procedure is qualified to use this device.
- This device must not be in contact with any electric instruments.

STORAGE

It is advisable to store the device in its original packaging at: Room temperature.

Use before the expiry date specified on the packaging.

POTENTIAL COMPLICATIONS

Endoscopic urological stone elimination procedures should only be performed by physicians familiar with the inherent potential complications of this type of procedures. Complications may occur during or after the procedure, and may include the following:

- Edema.
- Dissection.
- Urosepsis.
- Perforation of the ureter.
- Iatrogenic infections.

PROCEDURE

PREPARATION:

- Make sure the product packaging is not damaged. Do not use the product to perform the procedure if the packaging is damaged.
- Before the procedure, remove the product from its packaging in a suitable area, taking the necessary precautions for handling sterile products.
- Attach the Luer-lock valve to the optical system.
- Introduce the optical system into the urinary tract.
- Maneuver the closed retrieval basket through the optical system towards the stone, so the end basket is positioned in front of the stone.
- Trigger the handle control in order to open the retrieval basket.
- Maneuver the device to trap the stone in the retrieval basket.
- When the stone is inside the basket, trigger the handle control to close the retrieval basket in order to trap the stone.
- Remove the trapped stones.

NOTE: After removing the handle, it is safe to remove the sheath through the optical system, while the retrieval basket remains in place. See Handle removal.

HANDLE REMOVAL:

Press the button on the back of the device and remove the control from the handle. Using the distal end of the handle, turn the control screw in the direction opposite to that indicated by the arrows. Remove the sheath.

SYMBOLS USED IN LABELS

	PRODUCT CODE
	LOT NUMBER
	SERIAL NUMBER
	CAUTION
	DO NOT RESTERILIZE
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	READ THE INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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Fabricado por: PROMEDON S.A.
PROMEDON S.A.