# URODRAIN

Guidewire / Cable guía / Fio guia

M-40-15-1692 / PSD-00000043 (04) / SEP-2022

**User Manual**  
ENGLISH

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**Promedon**

## USER MANUAL

### ATTENTION!

Sterile product for single use (sterilized with ethylene oxide). Do not reuse or resterilize.

Do not use the device if packaging is open or damaged.

The device should be stored in its original packaging at ROOM

TEMPERATURE.

USE BEFORE THE EXPIRATION DATE SPECIFIED ON THE PACKAGING.

This product should be used only by physicians who are familiar with its application.

### DESCRIPTION

The diameter and length dimensions of the guidewire are specified on the label attached to the product packaging.

The guidewire type and material are specified on the label attached to the product packaging.

### PACKAGING CONTENT

The packaging has a guidewire whose dimensions are specified on the label attached to the product packaging.

Hydrophilic guidewires come in a protective plastic tube inside the product packaging.

### INSTRUCTIONS FOR USE

These guidewires are designed to facilitate placement, procedures and exchange of devices during endourological procedures.

### CONTRAINDICATIONS

Do NOT use this guidewire in the vascular system.

Check specific additional contraindications on the label of the device to be used with this guidewire.

### PRECAUTIONS

Check instructions for use, contraindications and possible complications in the instructions provided with the devices to be used with this guidewire.

This device must be used only by physicians trained in endourological procedures.

Check the guidewire carefully for kinks, L-shapes or other defects before use and, if possible, during the procedure.

Do not use damaged guidewires as trauma may be caused.

Before use, make sure the guidewire diameter is compatible with the device to be used.

The free movement of the guidewire inside the device to be used is an important feature of the directional guidewire system as the user is provided with valuable tactile information. Check the system for resistance before use.

Do not twist the guidewire without checking the movement of its tip; otherwise, trauma may be caused.

Twisting a guidewire offering resistance may damage the device.

Always advance or remove the guidewire slowly. Do not push or twist guidewires offering resistance. Resistance may be felt and/or noticed using an imaging system and observing if the tip bends.

Identify the cause of resistance and take action to eliminate it.

### PREPARATION

1. Follow the manufacturer's instructions to prepare the device to be used.

2. Wash the guidewire lumen before insertion.

3. For hydrophilically-coated guidewires, connect the syringe with saline solution to the Luer-lock connector and fill the protective jacket, moistening the full length of the guidewire. After checking that the full length of the guidewire has been moistened, remove it from its protective jacket for use.

### PROCEDURE

Check the instructions of the device to be used with the guidewire.

1. Insert the guidewire through the natural passage or the passage created according to the surgical procedure chosen.

2. Make sure the guidewire is in place using the imaging system.

3. Keep the guidewire in place during the procedure.

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**Fabricado por: PROMEDON S.A.**
PROMEDON S.A.