# CATÉTER URETERAL

M-70-15-3045 (02) / 23-Jun-2021

**Instructions for use**  
ENGLISH

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**Promedon** 

## INSTRUCTIONS FOR USE

### DESCRIPTION

The Promedon ureteral catheter is a single-use product used for drainage or retrograde pyelograms.

### PACKAGING CONTENTS

A radiopaque polyurethane catheter (see dimensions on the packaging label)

### INSTRUCTIONS

The Promedon ureteral catheter is prescribed for:

- Temporary urinary drainage
- Retrograde ureteropyelography
- Opaqing of the collecting system for percutaneous access
- Assistance in accessing the ureteral orifice
- Replacing guidewires
- Obtaining the cytology and brush biopsy of the upper tract
- Strengthening the proximal guidewire to enable distal handling

### CONTRAINDICATIONS

Stenoses preventing the placement of the catheter.

### WARNINGS

The catheter should not be handled with pointed, serrated or sharp objects, since any depression, deterioration, puncture or tear may cause subsequent complications. Fluff, fingerprints, powder or any other contaminated area may cause reactions to foreign bodies. Every precaution must be taken to avoid contamination. The bending or torsion of the catheter before or during placement can damage its integrity.

### PRECAUTIONS

The complete process must be performed under radiological control. Forcing the insertion of the catheter should be avoided.

Single-use sterile product (ethylene oxide-sterilized). It must not be reused or resterilized, since this could damage the device's performance and increase the risk of inadequate resterilization and cross contamination. Do not use the device if the packaging has been opened or damaged. The product is meant for use by a qualified physician who is familiar with ureteral endoscopic procedures.

### STORAGE

Store the device in its original packaging at ROOM TEMPERATURE.  
DO NOT USE AFTER THE EXPIRATION DATE ON THE PACKAGING.


### POTENTIAL COMPLICATIONS

- Ureter perforation

### PROCEDURE

1. Wash the ureteral catheter, the pusher, and the guidewire with saline solution.
2. Identify the flexible end of the guidewire.
3. Insert the flexible end of the guidewire through the work channel of the cystoscope to be used.
4. Locate the ureter and introduce the guidewire until its proximal end reaches the renal pelvis.
5. Insert the ureteral catheter concentrically to the previously inserted guidewire.

SYMBOLS USED IN LABELS

	CATALOG NUMBER
	LOT NUMBER
	CAUTION
	SERIAL NUMBER
	EXPIRY DATE
	DO NOT REUSE
	DATE OF MANUFACTURE
	DO NOT RESTERILIZE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	MANUFACTURER
	READ THE INSTRUCTIONS FOR USE
	DO NOT USE IF PACKAGING IS DAMAGED
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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**Fabricado por: PROMEDON S.A**
PROMEDON S.A.