

vontris[®]
VUR treatment

A DEFINITIVE SOLUTION FOR VUR



RELIABLE RESULTS,
EVEN IN HIGH-GRADE REFLUX

Promedon
— People + Innovation

Endoscopic treatment of vesicoureteral reflux

Preferred for the management of Vesicoureteral Reflux (VUR) in children, endoscopic treatment has numerous advantages over long-term antibiotic prophylaxis and surgical procedures.

Through a subureteral injection of tissue bulking substance, endoscopic treatment of VUR offers an effective and minimally invasive alternative to surgery when antibiotic prophylaxis fails to produce the desired results.



The benefits of endoscopic treatment


- Minimally invasive procedure.
- Treatment performed on an outpatient basis.
- Immediate results.
- Higher efficacy rate than antibiotic prophylaxis.
- No surgery-related risks.
- Unlike the surgical alternative, hospital visits and costs are significantly reduced.



Safe and Effective Treatment

Vantris is a biocompatible, non-absorbable, synthetic tissue-bulking agent designed for the treatment of Vesicoureteral Reflux (VUR) in children.

As a hydrogel comprised of stable, highly malleable macroparticles, Vantris is resistant to particle migration.



Why Choose Vantris?



- High success rate --96,9% to 97,77%-- with an average injection volume in a range from 0,4 ml to 0,7 ml. Mean follow-up of up to 54 months.⁽¹⁻²⁻³⁾
- Long-term bulking effect.
- Vantris has demonstrated to be effective even in high-grade reflux and also in complex cases.⁽³⁻⁴⁻⁵⁻⁶⁾
- According to the reports of the studies performed, absence of migration is observed mainly due to the size of its macroparticles, an average of 300 microns.
- Vantris is compound of biocompatible material that is non-immunogenic and non-antigenic. Additionally, its non-animal origin greatly reduces the risk of an allergic reaction⁽⁷⁾.
- Precise injection due to the high-fluidity hydrogel that consists of amorphous and flexible macroparticles that can be extruded with 22 and 23-gauge needles.

According to the description of an ideal tissue bulking agent ⁽⁸⁾, Vantris embodies the necessary characteristics to achieve a safe, effective and long-lasting result.


Vantris: substance properties

Vantris is a substance consisting of polyacrylate/polyalcohol copolymer (PPC) particles immersed in a carrier that contains 40% glycerol.


Once implanted, the carrier is eliminated unmetabolized by the reticuloendothelial system and excreted through the kidneys. However, the Vantris particles remain to ensure long-lasting bulking.



Injection day.



12 Weeks.



2 Years.

- Structural water
- Glycerol
- Tissue
- PPC

The safe and innovative technology of Vantris

- Vantris has been developed as a non-absorbable, synthetic, injectable material, with high stability and durability, aiming at long-term effectiveness.
- Vantris hydrogel particles are comprised of a polyacrylate/polyalcohol copolymer.



- In recent years, the components of Vantris have been successfully and widely utilized throughout the medical and biotechnology industries for the following:

- Intraocular lenses
- Artificial organs
- Injectable material for plastic and reconstructive surgery
- Drug delivery systems
- Bulking for embolization of hypervascularized tumors

- Biocompatibility: The material has been successfully tested in cytotoxicity, sensitization, irritation or intracutaneous reactivity, acute systemic toxicity, subchronic and subacute toxicity, genotoxicity, implantation, chronic toxicity, carcinogenicity and migration tests.


- Vantris physicochemical properties make it a highly stable material that withstands thermal or pH changes that can arise in treated tissues.

- Vantris particles comply with the requirements of an **ideal biomaterial**⁽⁸⁾:

- Non-toxic, non-pyrogenic, non-hemolytic, non-inflammatory.
- Non-allergenic, non-carcinogenic, non-teratogenic, non-cytotoxic and painless for patients.
- Effective: functional, absence of migration, reliable, durable, and easy to implant material.
- Biocompatible.

- Biological tests have shown absence of migration to other organs, due to the size of its particles⁽⁹⁾; it produces no allergic reactions or chronic inflammatory processes (granuloma formation).

- Vantris has demonstrated, in a comparative evaluation of histopathological changes through time with two commonly used substances, to be the bulking agent that generates the least tissue reaction and inflammatory infiltration, minimum fibrotic tissue and also proved that it does not generate any calcifications.⁽¹⁰⁾



Macroparticles

- When compressed, the malleable macroparticles of Vantris change in shape to enable extrusion with **22 and 23-gauge** needles.
- When these particles are implanted in the ureterovesical junction, the material acts as an enlarger, increasing the volume of the area and correcting the anatomy of the meatus and the distal ureter, preventing urine from returning to the ureter after being stored in the bladder.
- The particles increase the tissue volume, generating a minimum fibrotic growth around them, 70 microns thick.




Superior Performance

Vantris macroparticles have an average size of 300 μm , which is a key factor in the absence of migration reported in literature. (Fig. 1).

The testing of Vantris in animals has proven the absence of particle migration.⁽⁹⁾

Fig. 1: Comparison of injectable material particles size.




The risk of migration is directly correlated with the size of the particles. (Fig. 1).

An effective, simple, and safe procedure

As a minimally invasive procedure, an endoscopic injection of Vantris requires no hospitalization and is usually free of any short-term and long-term complications.

In most cases, this outpatient procedure results in a quick recovery from the anesthesia as well as an immediate correction of the VUR. This in turn allows patients to return to their normal daily activities after a short period of time.

Vantris provides a high level of reflux resolution after single endoscopic injection. Due to the substance features, Vantris has demonstrated high performance with very low VUR recurrence, unlike other biodegradable bulking agents which show up to 26% reflux recurrence after one year.⁽¹¹⁾



Step-by-step STING procedure:

In most cases, Vantris is injected under general anaesthesia.


In order to perform the procedure, the following components are necessary:

- 1 syringe. The quantity to be used will depend on each particular case.
- 1 flexible injection needle or alternatively 1 metal semi-rigid injection needle. It is recommended to use the injection needles developed by Promedon, with the following features:
 - 3.7 Fr or 5 Fr Flexible Injection Needle, 23 G Bevel tip.
 - 3.6 Fr Metal Semi-rigid Injection Needle with 22 G Bevel tip or Concave side opening tip.
- 1 cystoscope with a straight working channel of 4 Fr or more, according to the chosen needle.

1. The free flow through the injection needle is verified using saline solution.
2. The syringe is connected to the injection needle passing material through it until the substance appears at the needle tip. Insert the injection needle into the cystoscope.
3. Prior to the injection, the meatus should be observed with different volumes of bladder filling in order to choose the ideal situation (Figure A).

Ensure that the tip of the needle is facing the ureteral side at a 6 o'clock position.

4. The submucosa of the bladder is punctured at the six o'clock position, 3 mm under the ureteral meatus, 4 to 5 mm depth (Figure B). The anatomy of the ureteral meatus defines the selection of the puncture site and number of punctures needed.
5. The material is then slowly injected until the ureteral wall is adequately modified (bulkiness). More punctures may be performed until the desired effect is reached (Figure C and D).
6. Once the injection is completed, the needle is kept in its position for 30 seconds. Then, the needle is removed.
7. Once the procedure has been completed, the bladder is emptied, the cystoscope removed, and an optional voiding cystoureterogram (VCUG) may be performed to verify reflux repair.




[Figure A] Ureteral meatus prior to injection.




[Figure B] Insertion site for needle.



[Figure C] Location of needle during injection.



[Figure D] Tissue growth in the meatus after a successful injection.



[Figure E] Final view of bulking in the ureteral meatus.



References:

1. Boris Chertin , Wael Abu Arafeh, Alexander Zeldin, Israel A. Ostrovsky, Andis Stanislav Kocherov. Department of Pediatric Urology, Shaare Zedek Medical Center, Faculty of Medical Science, Hebrew University, Jerusalem, Israel. Endoscopic correction of VUR using Vantris as a new non-biodegradable tissue augmenting substance. Three years of prospective follow-up. *Urology*, 2013 Feb 25 - doi: 10.1016/j.urology.2013.01.024.
2. Juan Pablo Corbetta*, Juan I. Bortagaray, Santiago Weller, Javier Ruiz, Carol Burek, Cristian Sager, Víctor Durán, Juan C. Lopez. Urology Department, Hospital de Pediatría Prof. Dr. Juan P. Garrahan, Argentina.
The use of polyacrylate-polyalcohol copolymer hydrogel in the endoscopic treatment of primary vesicoureteral reflux in children. *Journal of Pediatric Surgery*. DOI: <http://dx.doi.org/10.1016/j.jpedsurg.2014.07.001>
3. Stanislav Kocherov, Ibrahim Ulman, Sergey Nikolaev, Juan Pablo Corbetta, Rudin YE Zafer Dokumcu, Ali Avanoğlu, Ludmila Menovshchikova, Semen Kovarskiy, Tatiana Skliarova, Santiago Weller, Juan I Bortagaray, Juan C Lopez, Víctor Durán, Carol Burek, Cristian Sager, Maruhnenko DV, Garmanova TN, Aliev DK, Wael Abu Arafeh, Boris Chertin. Multicenter survey of endoscopic treatment of vesicoureteral reflux using polyacrylate-polyalcohol bulking copolymer (Vantris). *Urology*. 2014 Sep;84(3):689-93. doi: 10.1016/j.urology.2014.04.033.
4. Boris Chertin, Wael Abu Arafeh, Stanislav Kocherov. Department of Pediatric Urology, Shaare Zedek Medical Center, Faculty of Medical Science, Hebrew University, Jerusalem, Israel. Endoscopic correction of complex cases of vesicoureteral reflux utilizing Vantris as a new non-biodegradable tissue-augmenting substance. *Pediatr Surg Int*. 2014 Jan 22. [Epub ahead of print]. DOI 10.1007/s00383-014-3468-z
5. Maria Ormaeche, Juan Moldes, Eduardo Ruiz, Cesar Benmaor, Ricardo Soria, Andres Villegas, Roberto Vagni, Francisco Debadiaola. Hospital Italiano de Buenos Aires, Buenos Aires. Argentina. Manual Injectable Treatment for High-Grade Reflux with a Permanent Bulking Substance (Vantris). 1st World Congress of Pediatric Urology. San Francisco, USA. May 2010. Abstract# 33.
6. Ali Tekin, Ismail Yagmur, Ulkum Zafer Dokumcu, Raziye Ergun, Ali Avanoğlu and Ibrahim Ulman. Ege University Faculty of Medicine, Division of Pediatric Urology, Department of Pediatric Surgery, Izmir, Turkey. Endoscopic treatment of Vesicoureteral Reflux using two different bulking agents: A novel material polyacrylate polyalcohol copolymer vs. Dextranomer Hyaluronic Acid. Abstract S2-4. ESPU Congress 2014, Innsbruck, Austria.
7. Data on file. Biocompatibility Tests. Promedon, 2004/2005.
8. Dmochowski RR, Appell RA. Injectable agents in the treatment of stress urinary incontinence in women: where are we now?. *Urology*. 2000 Dec 4;56(6 Suppl 1):32-40.
9. Ormaechea M, Paladini M, Pisano R, Scagliotti M, Sambuelli R, Lopez S, Guidi A, Muñoz J, Rossetti V, Carnerero M, Beltramo D, Alasino R, Bianco I, Griguol O, Valladares D, De Badiola F. Universidad Católica de Córdoba, School of Medicine, Córdoba, Argentina. Vantris, a biocompatible, synthetic, non-biodegradable, easy-to-inject bulking substance. Evaluation of local tissular reaction, localized migration and long-distance migration. *Arch Esp Urol*. 2008 Mar;61(2):263-8.
10. De Badiola F, Villegas Scivetti A., Soria R., Vagni R., Centurion J., Ormaechea M., Moldes J., Ruiz E. Hospital Italiano de Buenos Aires, Buenos Aires – Argentina. Histopathological changes after Polyacrylate Polyalcohol Copolymer (Vantris) injection. CAU-SIUP Meeting Chile 2010. Santiago, Chile. Abstract# PPC6. [Spanish].
11. Lee EK, Gatti JM, DeMarco RT, Murphy JP. Long-term follow-up of dextranomer/hyaluronic acid injection for vesicoureteral reflux: late failure warrants continued followup. *J Urol* 2009; 181:1869e74.

Ordering information

Product Code: BAR 1J

1 ml Vantris syringe – Ref: BARI-1J

