



Testicular Prosthesis

NATURAL AESTHETIC

The availability of different sizes makes it possible to restore the scrotal aesthetic to the patient's individual anatomical circumstances.

SAFE

The soft-solid silicone design prevents failure by leakage, loss of shell integrity or shrinking over time.




Testicular Prosthesis

N&S is a soft-solid implant made of low hardness silicone elastomer, to resemble natural aesthetic of the testicle. It is designed to replace the testicle in shape and size in one or both sides to restore the normal appearance of the scrotum and improve the Quality of Life of the patient.

N&S is indicated in patients who have suffered from testicular agenesis or for those who have had or need to have a testicle removed through surgery due to different pathologies.

ORDERING INFORMATION



Ref. Number	T-Small	T-Medium	T-Large	T-XLarge
A	31 mm	37 mm	42 mm	47 mm
B	23 mm	28 mm	32 mm	34 mm
Volume	9 ml	15 ml	23 ml	28 ml

CONTRAINDICATIONS

N&S must not be implanted in:

- Patients who present an infection in any part of the body, particularly if it is genital or in the urinary tract.
- Patients with neoplasia or transscrotal infection.

CE 0197 PSD-0000053 (07) / JULY-2024